Baerveldt™

Glaucoma Implant
Implant De Drainage Du Glaucome
Glaukom-Implantat
Glaucom Implantaat
Implante Para Glaucoma
Impianto Di Drenaggio Per Glaucoma
Glaukom Implantat
Glaukomimplantat
Implante Para Glaucoma
Glaukoomaistute
Εμφύτευμα γλαυκώματος

Cover front

title : Baerveldt
no. : 00084024
revisieno. : 04
date : 14 apr 2006
printcolor : black
text : black

checked

date

 FIELD BOXMORE HEALTHCARE PACKAGING	Your item code:	00084024/04	Colours used
	Item code Mareen:	i81945.eps	● ZWART 00013
	Cutting profile No.:	75851	● KAP
	Dimensions:	150X100 mm	
	Date:	14 apr 2006	
	Proof by:	Burick Christophe	
	Proof No.:	1	
 Chesapeake Company			

Caution: Federal (USA) law restricts this device to sale by or on the order of a physician.

GLAUCOMA IMPLANT/DEVICE DESCRIPTION: The Baerveldt™ Glaucoma Implant and the Baerveldt™ Pars Plana Glaucoma Implant are non-valved, artificial drainage devices that shunt aqueous via a tube to an episcleral plate centered over the equatorial region of the globe. They are constructed entirely of medical grade silicone with an episcleral plate with different surface areas of between 200 mm² and 500 mm². The plate of a Baerveldt™ implant is molded to conform to the curvature of the globe with a 13 mm convex radius. The silicone plate is barium impregnated to increase ultrasound resolution and identification with the CT Scan, MRI, and plain skull films. This also allows for the construction of a thin, durable, flexible implant. The proximal portion of the plate has flange with two large fixation holes for scleral attachment. The external diameter of the silicone tube is 0.63 mm and the internal diameter is 0.30 mm.

The Baerveldt™ Implant contains a silicone tube (32 mm in length) attached to the anterior surface of the plate. This tube is for insertion into the anterior chamber.

The Baerveldt™ Pars Plana Implant contains a silicone tube (7 mm in length) connected to a "elbow" attached to the anterior surface of the plate,

This elbow is for insertion into the posterior chamber.

Model	Surface Area	Plate Length	Tube Length	Insertion	Style
BG101-350	350 mm ²	32 mm	32 mm	Anterior Chamber	Straight Tube
BG102-350	350 mm ²	32 mm	7 mm	Posterior Chamber	90° Elbow Pars Plana
BG103-250	250 mm ²	22 mm	32 mm	Anterior Chamber	Straight Tube

INDICATIONS: For use in patients (with prior vitrectomy for Pars Plana) with medically uncontrollable glaucoma and poor surgical prognosis, such as, but not limited to: neovascular glaucoma, aphakic/pseudophakic glaucomas, patients who have failed conventional surgery, congenital glaucomas and secondary glaucomas due to uveitis, epithelial downgrowth, etc.

WARNINGS: Do not use the device if sterile package integrity has been compromised. Do not resterilize the implant by any method. Do not reuse the implant. Do not store at temperatures above 45°C (113°F).

CONTRAINDICATIONS: Bacterial conjunctivitis, bacterial corneal ulcers, endophthalmitis, orbital cellulitis, bacteremia or septicemia, active scleritis and/or no light perception.

COMPLICATIONS/ADVERSE EVENTS: The complications during and after surgery include, but are not limited to: choroidal hemorrhage, hyphema, serous choroidal effusion, hypotony, flat anterior chamber, phthisis bulbi, retinal detachment, endophthalmitis, tube erosion, tube touch to cornea, tube block by iris or vitreous, bullous keratopathy, uveitis and diplopia.

BACKGROUND PANEL
PRINTS 20% BLACK
(TYPICAL)

title : Baerveldt
no. : 00084024
revisieno. : 04
date : 14 apr 2006
printcolor : black
text : black

checked

date

	Your item code:	00084024/04	Colours used
	Item code Mareen:	i81945.eps	● ZWART 00013
	Cutting profile No.:	75851	● KAP
	Dimensions:	150X100 mm	
	Date:	14 apr 2006	
	Proof by:	Burick Christophe	
	Proof No.:	1	
			

DIRECTIONS FOR USE: 1. Examine the implant prior to implantation for proper size, model, and expiration date. 2. Open the pouch and remove the implant by cutting the suture in a sterile environment. 3. The implant can become electrostatically charged upon opening the package. Examine the implant to ensure that particles are not present on it. Rinse implant in sterile saline, if required.

IMPLANTATION: A high level of surgical skill is required for the implantation of artificial drainage devices. Correct patient selection and meticulous surgery and postoperative care are required. The implanting surgeon should have observed and/or assisted in numerous artificial glaucoma drainage device implantations prior to implanting a Baerveldt™ or a Baerveldt™ Pars Plana Glaucoma Implant. For both types of implants the tube should be ligated with a biodegradable ligature or obstructed with a temporary plug.

HOW SUPPLIED/ EXPIRATION DATE: Each Baerveldt™ and Baerveldt™ Pars Plana Glaucoma Implant is supplied sterile, in the dry form, in a container on a 26 mm diameter hemisphere with the implant attached to it by a 9-0 suture to ensure that the implant maintains its contours. The container is sealed within a double sterile pouch pack. The implants are gamma sterilized. Sterility is guaranteed unless the pouch is damaged or otherwise compromised. Expiration date is indicated on the outside of the box. The implant should not be used after the expiration date indicated.

RETURN OF DAMAGED/UNSTERILE PRODUCT: Contact Advanced Medical Optics, Inc. (AMO) for product return policy. Return the implant with proper identification and reason for return. Label return as a biohazard.

PATIENT INFORMATION: It is recommended that each patient receive information regarding Glaucoma implants prior to the decision to implant.

PATIENT CARD: An ID card is supplied in the implant packaging. This card should be given to the patient with instructions to keep as a permanent record of the implant and to show the card to any eye care practitioner seen in the future.

REPORTING: Adverse reactions and/or potentially sight-threatening complications that may be reasonably regarded as product related and that were not previously expected in nature, severity or incidence must be reported to AMO. This information is requested from all surgeons in order to document potential long-term effects of glaucoma implants. Potential problems may be related to a specific lot of product or may be indicative of long-term problems associated with these types of products. Report any product related event to AMO (800) 366-6554 U.S.A. Outside of the U.S. contact your local company representative.

DISCLAIMER OF LIABILITY: AMO accepts no liability for the choice of method or technique to implant the product or for the choice of the product for a particular patient or patient's condition.

BAERVELDT and the AMO ADVANCED MEDICAL OPTICS Logo are trademarks of Advanced Medical Optics, Inc. Manufactured by AMO Groningen BV, 9728 NX Groningen, The Netherlands. Advanced Medical Optics, Santa Ana, CA 92705, USA Toll-free (800) 366-6554. The CE marked Glaucoma Implants comply with the European Council Directive 93/42/EEC of June 14, 1993.

 FIELD BOXMORE HEALTHCARE PACKAGING	Your item code:	00084024/04	Colours used
	Item code Mareen:	i81945.eps	● ZWART 00013
	Cutting profile No.:	75851	● KAP
	Dimensions:	150X100 mm	
	Date:	14 apr 2006	
	Proof by:	Burick Christophe	
	Proof No.:	1	
			

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